

COLOR-CODED PE-PERFUSION LINES

- ✓ color-coded lines reduce risk of confusion
- ✓ polyethylene
- ✓ lipid and alcohol resistant
- ✓ pressure proof up to 4,5 bar
- ✓ suitable for all standard infusion pumps (syringe pumps)

dimensions		product code	
length ¹	150 cm	green	0222310
inside-ø	1 mm	blue	0222311
outside-ø	2 mm	yellow	0222312
box	200 pcs	red	0222283

¹ Available in different lengths.



PERFUSION LINES



green

blue

yellow

red

! #&! %&() *+,-./:;<=>?@A# 71* #15 B0#%00DD:2 #2 %&8% H4%\$7(>=JK@LJLME(6 H4%& N|>=JK@LJLME4

PE-PRESSURE LINE

- ✓ male/female Luer lock connector

dimensions	
length ¹	15, 30 cm
inside-ø	1 mm
outside-ø	2 mm
box	200 pcs

product code	
15 cm (red)	025861
30 cm (blue)	025872
30 cm (red)	022285

INFUSION LINE KMT

- ✓ special length for more freedom of movement during bone marrow transplantation
- ✓ transparent
- ✓ suitable for all standard infusion pumps (syringe pumps)
- ✓ pressure proof up to 4,5 bar
- ✓ PVC-free
- ✓ lipid and alcohol resistant

dimensions	
length ¹	900 cm
inside-ø	1,5 mm
outside-ø	2,8 mm
box	50 pcs

product code	
04OR120900PF	
	cm

1 Other lengths available on demand



PRESSURE AND INFUSION LINES



Pressure line red



Infusion line KMT

Venous connection line with Needle Free 3-way stopcock

- ✓ PVC-free
- ✓ with blue 3-way stopcock and integrated Needle Free Port

Needle Free 3-way stopcock

- ✓ blue marked
- ✓ PVC-free
- ✓ with integrated Needle Free Port

Venous connection line


Inside-ø	3 mm
Outside-ø	4,1 mm
Length	15 cm
Product Code	02MD20615PU
box	50 pcs

Needle Free 3-way stopcock


Product code	02MM51013
box	50 pcs



Extensions



(1)



(2)

Venous connection line (1) and Needle Free 3-way stopcock (2)

**ARTERIAL EXTENSION LINE
WITH 3-WAY STOPCOCK**

- ✓ mini extension line with a glued red 3-way stopcock
- ✓ pressure-tested (up to 65 psi)

Dimensions	
length	15 cm
Inside-ø	1 mm
Outside-ø	2 mm
box	100 pcs

Product code
028360

**VENOUS EXTENSION LINE
WITH BLUE 3-WAY STOPCOCK**

- ✓ mini extension line with a glued blue 3-way stopcock
- ✓ pressure-tested (up to 65 psi)

Dimensions	
length	25 cm
Inside-ø	1 mm
Outside-ø	2 mm
box	100 pcs

Product code
028113

**RED CONNECTION LINE
WITH RED 3-WAY STOPCOCK**

- ✓ connection line with a glued red 3-way stopcock
- ✓ pressure-tested (up to 65 psi)

Dimensions	
length	15 cm
Inside-ø	1,5 mm
Outside-ø	3 mm
box	100 pcs

Product code
02MD40153A

**EXTENSIONS**

arterial extension



venous extension



connection line

**INFU-STAR
INFUSION-MANIFOLDS**

- ✓ transparent
- ✓ colour coded cocks
- ✓ fat- and alcohol-resistant
- ✓ pressure-tested (up to 65 psi)
- ✓ PVC-free
- ✓ special lengths for increased mobility

Dimensions	
colour coded cocks	3, 5
extension line lengths	150, 250, 400, 900 cm
box	50 pcs

Product code	
3x manifold w/o connection line	028102F
5x manifold w/o connection line	028103F
3x manifold 250 cm	028116F
5x manifold 150 cm	02520612F
3x manifold 400 cm	02520614F
5x manifold 900 cm	025312
5x manifold with light protected line 900 cm	025312LI

MANIFOLDS-HOLDER

- ✓ secure fixation of manifolds on infusion stands
- ✓ suitable for all usual infusion stands
- ✓ suitable for steam sterilization
- ✓ fat- and alcohol-resistant

Product code
02AK6540

box	10 pcs
-----	--------




INFU-STAR



3x with / w/o connection line



5x w/o connection line



5x with connection line



holder

MiroMed® Pfeiderer Medical Produkte GmbH · Lindleystr. 12 · 60314 Frankfurt · Internet: www.miro-med.de · E-Mail: mail@miromed.de · Telefon (069) 489878-0 · Telefax (069) 489878-70

INFUSION EXTENSION LINE

- ✓ Luer-Lock hub

Dimensions	
length	75,100,140 cm
diameter	3,0 x 4,1 mm
box	100 pcs

Product code	
75 cm	02OR110100
100 cm	02OR112101
140 cm	02OR111100

PRESSURE CONNECTION LINE

- ✓ Luer-Lock male/male
- ✓ without stripes
- ✓ pressure-tested (up to 507 psi)

Dimensions	
length	100 cm
diameter	1,65 x 3,2 mm
box	100 pcs

Product code	
	02MD41002MM

PRESSURE CONNECTION LINE

- ✓ Luer-Lock male/female
- ✓ colour coded
- ✓ pressure tested (up to 507 psi)


Abmessungen	
length	150 cm
diameter	1,65 x 3,2 mm
box	50 pcs

Product code	
red (arterial)	02MD41503
blue (venous)	02MD41504

cm



INFUSION ASSECORIES



Infusion extension line

pressure connection line

pressure extension line blue

IRRIGATION SET

Closed system for irrigation and diagnostic radiology.

The set consists of:

- ✓ contrast bag 2 l, blue, with a large filler plug and connecting tube
- ✓ rectal tube with ballon blocking (maximal 90 ml air)
- ✓ syringe 50 ml with Luer Lock connector
- ✓ Latex-free

dimensions	
bag volume	2000 ml
syringe volume	50 ml
ballon capacity	90 ml
box	20 pcs

product code

02050022665

**IRRIGATION SET**

Irrigation set

TRANSFUSION SET P

Non-vented blood transfusion set

- ✓ transparent filter body
- ✓ type «P» for pressure transfusion
- ✓ for 1 to 3 blood units
- ✓ 40 µm filter membrane

dimension

particle filter	40 µm
filter area	32 cm ²
tubing length	150 cm
box	100 pcs

product code

2OR01440 0

HAEMOMED
TRANSFUSION SET G

Non-vented blood transfusion set

- ✓ typ «G» for gravity infusion
- ✓ 200 µm filter membrane

dimension

particle filter	200 µm
tubing length	150 cm
box	100 pcs

product code

02S60011



TRANSFUSION SET



Transfusion set