

MINI-CONNECTING TUBE FOR PEDIATRICS

- ✓ remaining volume less than 0,4 ml
- ✓ polyurethan
- ✓ lipid resistant Luer lock connector

dimensions	
length	25 cm
diameter	1,0 x 2,0 mm
box	100 pcs

product code	
male Luer lock connector	
light protected	058024LI
transparent	058024
male Luer Slip connector	
light protected	058026LI
transparent	058024

MINI-CONNECTING TUBE FOR PEDIATRICS

Characteristics similar to product above, additionally

- ✓ light protected

dimensions	
length	30 cm
box	100 pcs

product code
05MM31003LI

MINI-CONNECTING TUBE WITH INTEGRATED LIQUID FILTER FOR PEDIATRICS

- ✓ liquid filter 0,2 and 1,2 µm
- ✓ light protected line
- ✓ self-venting

dimensions		product code
length	150 cm	
filter membrane	0,2 µm	
diameter	1,0 x 2,0 mm	05MM31006LI
length	150 cm	
filter membrane	1,2 µm	
diameter	1,0 x 2,0 mm	05MM31005LI
box	50 pcs	

All products and packings are Latex- and PVC-free.



MINI-CONNECTING TUBES



transparent



light protected



with filter

! # \$ % & ' () * + , - . / : ; < = > ? @ [\] ^ _ ` { | } ~ ¡ ¢ £ ¤ ¥ ¦ § ¨ © ª « ¬ ® ¯ ° ± ² ³ ´ µ ¶ · ¸ ¹ º » ¼ ½ ¾ ¿ À Á Â Ã Ä Å Æ Ç È É Ê Ë Ì Í Î Ï Ñ Ò Ó Ô Õ Ö × Ø Ù Ú Û Ü Ý Þ ß à á â ã

MIRO-PED MICRO-CONNECTING TUBE FOR PEDIATRICS₁

- ✓ low remaining volume 0,2 ml/m
- ✓ lipid resistant luer lock connector

dimensions	
length	30-150 cm
diameter	0,5 x 2,0 mm
filling volume	0,2 ml/m
box	150 pcs

product code	
30 cm	058027
40 cm	058025
60 cm	058028
100 cm	058029
150 cm	058030

MIRO-PED MINI-Y-CONNECTOR₁

- ✓ light protected
- ✓ integrated check valves
- ✓ lipid resistant luer lock connector

ordering information		product code
box	50 pcs	05MM31001LI

MIRO-PED MINI-TRIPLE-CONNECTOR₁ similar to product 05MM31001LI, has additional feature

- ✓ needle free adapter

ordering information		product code
box	50 pcs	05MM31002LI

MINI-CONNECTING TUBE FOR PEDIATRICS

- ✓ with clamp
- ✓ low remaining volume 0,4 ml/m
- ✓ polyurethane
- ✓ lipid resistant luer lock connector

dimensions		product code
length	30 cm	05OR126300
diameter	1,0 x 2,0 mm	
box	100 pcs	

¹ All products and packings are
Latex- and PVC-free



MIRO-PED



Micro-connecting
tube



Mini-Y-Connector



Mini-Triple
Connector



Mini-Connecting
tube with clamp

INFUSION MANIFOLD FOR PEDIATRICS

- ✓ transparent body
- ✓ colored stop-cocks
- ✓ lipid and alcohol resistant
- ✓ pressure proof up to 4,5 bar
- ✓ PVC-free
- ✓ special length for more freedom of movement of the patient during therapy
- ✓ 0,2 µm transparent filter, with 96 hours of service time and automatic air-venting
- ✓ blue 1,2 µm particle filter with automatic air-venting
- ✓ light protected lines

dimensions	
colored stop-cock	3 pcs
lines length	100/150cm
box	25 pcs

product code

055476F



PEDIATRICS INFUSION SET



manifold with filters and connection lines

! #5! %8 () *+789##(%8, - €) #8 / 0/12 34 (5 67890##(-< (5 => ?; @/## 7/ * . #1(5 67890##(DD:2 "52 %8&8%6 H48\$87(!>=JK@LJLIME(6 H48% N|>=JK@LJLIME4

INFUSION SET MODEL SZT

Model SZT is suitable for drug therapy during stem cell transplantation.

- ✓ transparent body
- ✓ color coded stop-cocks
- ✓ lipid and alcohol resistant
- ✓ pressure proof up to 4,5 bar
- ✓ PVC-free
- ✓ 600 cm transparent spiral connection line
- ✓ three-way stop-cock (blue)

dimensions	
colored stop-cock	5 pcs
diameter	1,5 x 2,8 mm
lines length ¹	600 cm
box	50 pcs

product code

05MM41229

¹ Special length lines and light protected lines are available upon request.



INFUSION SET MODEL SZT

NEW PRODUCT



manifold with 5 stop-cocks and spiral connection line