

MiroMed® Pfeiderer Medical Produkte GmbH · Lindleystr. 12 · 60314 Frankfurt · Internet: www.miro-med.de · E-Mail: mail@miro-med.de · Telefon (069) 489878-0 · Telefax (069) 489878-70

FOLEY SILICONE BALLOON CATHETER WITH NELATON TIP

- ✓ best patient comfort
- ✓ reduced incrustation
- ✓ transparent
- ✓ x-ray contrast stripes
- ✓ single sterile packed catheter stoppers come with each catheter

Dimensions	
length	40 cm
diameter	12–26 CH
balloon capacity	5 ml, 30 ml
box	10 pcs

product code	
	08501205
	08501430
	08501630
	08501830
	08502030
	08502230
	08502430
	08502630

CH ml

FOLEY SILICONE BALLOON CATHETER WITH NELATON TIP FOR KIDS

- ✓ best patient comfort
- ✓ reduced incrustation
- ✓ transparent
- ✓ x-ray contrast stripes
- ✓ single sterile packed catheter stoppers come with each catheter

Dimensions	
length	33 cm
diameter	6–10 CH
balloon capacity	1 ml, 3 ml
box	10 pcs

product code	
	08500601
	08500803
	08501003

CH ml

FOLEY SILICONE BALLOON CATHETER WITH NELATON TIP, INTEGRAL

- ✓ best patient comfort
- ✓ reduced incrustation
- ✓ transparent
- ✓ x-ray contrast stripes
- ✓ single sterile packed catheter stoppers come with each catheter

1 integrated balloon type

Dimensions	
length	40 cm
diameter	12–24 CH
balloon capacity	5 ml, 10 ml
box	10 pcs

product code	
	08581205
	08581405
	08581610
	08581810
	08582010
	08582210
	08582410

CH ml



SILICONE BALLOON CATHETER



Nelaton Tip



Nelaton Tip, for kids



Nelaton Tip (integral)

FOLEY SILICONE BALLOON CATHETER WITH TIEMAN TIP

- ✓ best patient comfort
- ✓ reduced incrustation
- ✓ transparent
- ✓ x-ray contrast stripes
- ✓ single sterile packed catheter stoppers come with each catheter

Dimensions	product code
length 40 cm	08511003
diameter 10–20 CH	08511205
balloon capacity 3 ml, 5 ml, 30 ml	08511430
	08511630
box 10 pcs	08511830
	08512030
	CH ml

FOLEY SILICONE COATED LATEX CATHETER WITH NELATON TIP

- ✓ best patient comfort
- ✓ blue coloured
- ✓ best possible price performance ratio
- ✓ single sterile packed catheter stoppers come with each catheter

Dimensions	product code
length 40 cm	08521230
diameter 12–24 CH	08521430
balloon capacity 30 ml	08521630
box 10 pcs	08521830
	08522030
	08522230
	08522430
	CH ml

1 Don't use it if you are allergic to latex!

**BALLOON CATHETER**

Tieman Tip



Nelaton Tip

MiroMed® Pfeiderer Medical Produkte GmbH · Lindleystr. 12 · 60314 Frankfurt · Internet: www.miro-med.de · E-Mail: mail@miro-med.de · Telefon (069) 489878-0 · Telefax (069) 489878-70

3-WAY FOLEY SILICONE IRRIGATION CATHETER WITH DUFOUR TIP

- ✓ high balloon capacity
- ✓ transparent
- ✓ greater wall thickness
- ✓ catheter stoppers are individually sterile packed and come with each catheter

Dimensions	
length	40 cm
diameter	18–24 CH
balloon capacity	80 ml
box	10 pcs

product code	
	08531880
	08532080
	08532280
	08532480

CH ml

3-WAY FOLEY SILICONE IRRIGATION CATHETER WITH COUVELAIR TIP

- ✓ high balloon capacity
- ✓ transparent
- ✓ greater wall thickness
- ✓ catheter stoppers are individually sterile packed and come with each catheter

Dimensions	
length	40 cm
diameter	18–24 CH
balloon capacity	80 ml
box	10 pcs

product code	
	08541880
	08542080
	08542280
	08542480

CH ml

3-WAY FOLEY SILICONE IRRIGATION CATHETER WITH MERCIER TIP

- ✓ high balloon capacity
- ✓ transparent
- ✓ greater wall thickness
- ✓ catheter stoppers are individually sterile packed and come with each catheter

Dimensions	
length	40 cm
diameter	18–24 CH
balloon capacity	80 ml
box	10 pcs

product code	
	08551880
	08552080
	08552280
	08552480

CH ml



IRRIGATION CATHETERS



Dufour Tip



Couvelair Tip



Mercier Tip

Supracare Suprapubic Puncture Set Silicone Balloon Catheter¹

Suprapubic puncture set with balloon catheter.
The set consists of a silicone balloon catheter¹, a split trocar and a MiroValve® autom. catheter valve.

Silicone balloon catheter

- best patient acceptance
- reduced incrustation
- transparent
- open tip with side holes

Split trocar

- colour coded split trocar

MiroValve® autom. catheter valve

- pressure resistant against non controlled urine flow
- compatibel to all common bag connectors
- automatically open once connected to the bag
- enhanced patient comfort

1 integrated balloon type.



Supracare Technical Data



Length	40 cm
Diameter	12–16 FR
Balloon Capacity	5 ml (FR12+14), 10ml
Box	5 pcs.

Sizes

Product code

Split Trocar Length	12 cm	
Split Trocar-ø	5,6 mm	08980512
Catheter	12 FR	
Split Trocar Length	12 cm	
Split Trocar-ø	5,6 mm	08980514
Catheter	14 FR	
Split Trocar Length	15 cm	
Split Trocar-ø	6,3 mm	08980516
Catheter	16 FR	

Supracare Puncture Set



Silicone Balloon Catheter (1), Split Trocar (2), MiroValve®

MiroMed® Pfeiderer Medical Produkte GmbH · Lindleystr. 12 · 60314 Frankfurt · Internet: www.miro-med.de · E-Mail: mail@miro-med.de · Telefon (069) 489878-0 · Telefax (069) 489878-70

SUPRACARE CATHETER EXCHANGING SET WITH SILICONE SUPRAPUBIC / INDWELLING CATHETER AND EXCHANGE GUIDE WIRE

The set includes one silicone balloon catheter (integral) and one suprapubic exchange guidewire. The silicone balloon catheter has the following characteristics:

- ✓ best patient comfort
- ✓ reduced Incrustation transparent
- ✓ central & lateral hole

Dimensions	
length	40 cm
diameter	12–16 CH
balloon capacity	5 ml
box	5 pcs

product code	
12 CH	084312
14 CH	084314
16 CH	084316

CH

SUPRACARE, SUPRAPUBIC EXCHANGE GUIDE WIRE ONLY

The Suprapubic Exchange Guidewire has the following characteristics:

- ✓ straight flexible tip of 3 cm
- ✓ PTFE-coated


Dimensions	
length	80 cm
diameter	0,89 mm 0,035"
box	1 pcs

product code	
	089204

1 integrated balloon type



SUPRACARE EXCHANGE SET



Suprapubic Silicone balloon Cathter, Suprabubic Exchange Guidewire

SUPRACARE DILATOR EXCHANGE SET WITH SILICONE INDWELLING CATHETER, GUIDE WIRE AND SPLIT DILATOR

The set includes a silicone balloon catheter (integral), one guidewire and one dilator. The balloon catheter has the following characteristics:

- ✓ best patient comfort
- ✓ reduced incrustation
- ✓ transparent
- ✓ formed tip for dilation
- ✓ central lateral hole

Dimensions	
length	40 cm
diameter	12, 14 CH
balloon capacity	5 ml
box	10 pcs

GUIDE WIRE

The guide wire has the following characteristics:

- ✓ straight flexible tip of cm
- ✓ PTFE-coated

Dimensions	
length	80 cm
diameter	0,89 mm 0,035"

SPLIT DILATOR

The dilator has the following characteristics:


- ✓ peel away sheath PTFE coated
- ✓ atraumatic transition

Dimensions	product code
dilator-ø 14 CH catheter-ø 12 CH	084302
dilator-ø 16 CH catheter-ø 14 CH	084304

1 integral balloon type



DILATOR EXCHANGE SET



Suprapubic Silicone Balloon Catheter, Guidewire + Dilator

! #5! %8! () *+0288#/#! %8!, - +) #88 /0/612 34 (5 678&0899#!;< (5 => ?; @/## 7! * #! (5 670#7%000DD:2 #52 %8&8%/6 6F! - "C2 - "G 2 #52 %8&8%/6 6F#8\$7(!>=JK@LJLIME(6 H44%- N!>=JK@LJLIME4

**SUPRACARE UNIBAL SILICONE
INDWELLING CATHETER SINGLE,
STERILE**

The suprabubic silicone indwelling catheter features:

- ✓ no glue points between balloon and catheter shaft
- ✓ smooth distal tip
- ✓ best patient comfort
- ✓ reduced incrustation
- ✓ transparent
- ✓ central and lateral hole
- ✓ reinforced formed tip

dimensions	
length	40 cm
diameter	10–20 CH
balloon capacity	5, 10 ml
box	10 pcs

product code
08451005
08451205
08451405
08451610
08451810
08452010

CH ml

**SUPRACARE UNIBAL SILICONE
INDWELLING CATHETER SINGLE,
STERILE, SHORT**

The suprapubic silicone indwelling catheter features:


- ✓ best patient comfort
- ✓ reduced incrustation
- ✓ transparent
- ✓ central & lateral hole

dimensions	
length	28 cm
diameter	12–20 CH
balloon capacity	5, 10 ml
box	10 pcs

product code
08471205
08471405
08471610
08471810
08472010

CH ml

1 integrated balloon type

**SILICONE INDWELLING BALLOON CATHETER**

suprabubic silicone balloon catheter unibal & integral

MIROMARFLOW BLUE STENT
DOUBLE-J URETERAL STENT SET
TIEMANN ONE SIDE OPEN

The set consist of:

- ✓ ureteral stent made of aliphatic polyurethane
- ✓ pusher
- ✓ guidewire

technical data ureteral stent	
length	26-30 cm
diameter	6-8 CH
box	5 pcs

product code
08SS626G
08SS628G
08SS630G
08SS726G
08SS728G
08SS730G
08SS826G
08SS828G
08SS830G

CHcm

URETERAL STENT - SET
URETERAL STENT

- ✓ high quality aliphatic polyurethane
- ✓ perfect X-ray visibility
- ✓ low incrustation rate
- ✓ graduation in 5 cm steps
- ✓ high wearing comfort


URETERAL STENT - SET
GUIDEWIRE

- ✓ straight
- ✓ one side flexible tip
- ✓ PTFE-coated

technical data of guidewire	
length	125 cm
diameter	0,035" (by stents with 6 CH)
	0,038" (by stents with 7 and 8 CH)



DOUBLE-J URETERAL STENT SETS



Tiemann one side open stent (1) with pusher (2) and guidewire (3)

! # \$ % & ' () * + , - . / : ; < = > ? @ [\] ^ _ ` { | } ~ ¡ ¢ £ ¤ ¥ ¦ § ¨ © ª « ¬ ® ¯ ° ± ² ³ ´ µ ¶ · ¸ ¹ º » ¼ ½ ¾ ¿ À Á Â Ã Ä Å Æ Ç È É Ê Ë Ì Í Î Ï Ñ Ò Ó Ô Õ Ö × Ø Ù Ú Û Ü Ý Þ ß à á â ã

MIROMARFLOW BLUE STENT
MONO-J URETERAL STENT SET
TIEMANN, CLOSED

The set consist of:

- ✓ ureteral stent
- ✓ guidewire
- ✓ fixing clamp

technical data ureteral stent

length	70 cm
diameter	4,7-8 CH
box	5 pcs

product code

08MJC470GC
08MJC670GC
08MJC770GC
08MJC870GC

CHcm

URETERAL STENT - SET
URETERAL STENT

- ✓ high quality aliphatic polyurethane
- ✓ perfect X-ray visibility
- ✓ low incrustation rate
- ✓ graduation in 5 cm steps
- ✓ high wearing comfort

URETERAL STENT-SET
GUIDEWIRE


- ✓ straight
- ✓ one side flexible tip
- ✓ PTFE-coated

technical data of guidewire

length	70 cm
diameter	0,032" (by stents with 4,7 and 6 CH)
	0,038" (by stents with 7 and 8 CH)



MONO-J URETERAL STENT SETS



Tiemann, closed tip (1), fixing clamp (2) and guidewire (3)

MIROMARFLOW BLUE STENT
MONO-J URETERAL STENT SET
TIEMANN, BOTHSIDES OPEN

The set consist of:

- ✓ ureteral stent
- ✓ guidewire

technical data ureteral stent	
length	70 cm
diameter	4,7-8 CH
box	5 pcs

product code
08MJ0470G
08MJ0670G
08MJ0770G
08MJ0870G

CH cm

URETERAL STENT - SET
URETERAL STENT

- ✓ high quality aliphatic polyurethane
- ✓ perfect X-ray visibility
- ✓ low incrustation rate
- ✓ graduation in 5 cm steps
- ✓ high wearing comfort


URETERAL STENT - SET
GUIDEWIRE

- ✓ straight
- ✓ one side flexible tip
- ✓ PTFE-coated

technical data of guidewire	
length	70 cm
diameter	0,032" (by stents with 4,7 and 6 CH) 0,038" (by stents with 7 and 8 CH)



MONO-J URETERAL STENT - SET



Tiemann, open tip (1), guidewire (2)

! # \$ % & ' () * + , - . / : ; < = > ? @ [\] ^ _ ` { | } ~ ¡ ¢ £ ¤ ¥ ¦ § ¨ © ª « ¬ ® ¯ ° ± ² ³ ´ µ ¶ · ¸ ¹ º » ¼ ½ ¾ ¿ À Á Â Ã Ä Å Æ Ç È É Ê Ë Ì Í Î Ï Ñ Ò Ó Ô Õ Ö × Ø Ù Ú Û Ü Ý Þ ß à á â ã

NEPHROSTOMY EXCHANGE SET WITH BALLOON CATHETER (INTEGRAL)

The set includes a silicone balloon catheter (intergral), one guide wire (Lunderquist) and 3 dilators.

Dimensions	
length	40 cm
diameter	12, 14 CH
balloon capacity	5 ml
box	10 pcs

- ✓ best patient comfort
- ✓ reduced incrustation
- ✓ transparent
- ✓ x-ray contrast stripes
- ✓ central opening
- ✓ two lateral holes before the balloon
- ✓ one hole behind the balloon

GUIDE WIRE

Dimensions	
length	90 cm
diameter	0,89 mm 0,035"

- ✓ straight flexible tip of 3 cm
- ✓ PTFE-coated



DILATATOREN

	diameter	product code
dilator-ø	8, 10, 12 CH	Set 1
splitdilator-ø	14 CH	
catheter-ø	12 CH	081412
dilator-ø	10, 12, 14 CH	Set 2
splitdilator-ø	16 CH	
catheter-ø	14 CH	081414

1 integrated balloon type



NEPHROSTOMY EXCHANGE SET



Silicone balloon catheter

Guidewire

Dilator Set

! #& %& () *+0123456789:;<=>?@A#7/*.#15B0#%&00DD:2 "#2 %&8% H#%\$7(!>=JK@LJLIME(6 H#%& N|>=JK@LJLIMEA

PERCUTANEOUS NEPHROSTOMY SET PIGTAIL TYPE

This catheter is introduced according to Seldinger technique: coaxial or parallel. The set consist of:


- ✓ Pigtail-drainage-catheter, 30 cm length, lateral hole
- ✓ Initial-punction-needle, radiopaque, 18 G (1,27 mm), three parts
- ✓ PTFE Dilatators, 20 cm length
- ✓ PTFE-coated guide wire, 0,038", 3 mm J-Tip, 80 cm long
- ✓ Drainage tube (male Luer Lock connector), 20 cm length
- ✓ scalpel

dimensions		product code	
length	30 cm		08PDSP7
diameter	7-10 CH		08PDSP8
ø guidewire	0,038"		08PDSP9
box	5 pcs		08PDSP10

CH



NEPHROSTOMY SET



Nephrostomy set pigtail type

! #5! %8! () *+789##(%8!, - +) #8& /0/12 34 (5 67&8999#(-< (5 => ?; @/## 7/ * . #1) 5 B0/#%00DD:2 "#2 %&&8% H48\$87 (>=JK@LJLIME(6 H48%- N|>=JK@LJLIME4

PERCUTANEOUS NEPHROSTOMY
CATHETER PIGTAIL TYPE

- ✓ pigtail-drainage-catheter
- ✓ length 30 cm
- ✓ lateral opening

dimensions	
length	30 cm
diameter	7-10 CH
box	5 pcs

product code	
	08PCNP7
	08PCNP8
	08PCNP9
	08PCNP10

CH



NEPHROSTOMY CATHETER



catheter pigtail type

! #& % () * + , - . / : ; < = > ? @ [\] ^ _ ` { | } ~ ¡ ¢ £ ¤ ¥ ¦ § ¨ © ª « ¬ ® ¯ ° ± ² ³ ´ µ ¶ · ¸ ¹ º » ¼ ½ ¾ ¿ À Á Â Ã Ä Å Æ Ç È É Ê Ë Ì Í Î Ï Ñ Ò Ó Ô Õ Ö × Ø Ù Ú Û Ü Ý Þ ß à á â ã

RENODRAIN PERCUTANEOUS
NEPHROSTOMY SET
ACC. DR. SCHÜLLER

The set includes:

- ✓ Nephrostomy catheter, diameter 8 CH, 25 cm
- ✓ 3-part needle
- ✓ obturator
- ✓ universal adapter
- ✓ Guidewire (Lunderquist), J-tip, PTFE-coated

Dimensions	
length	25 cm
diameter	8 CH
box	10 pcs


product code
08RE400870

CH

! "#\$ %& '()*+,-./:;<=>?@A[B CDEFGHIJKL MNOPQRSTUVWXYZ[\]^_`{|}~!@#\$%^&*()
Pfliderer Medical Produkte GmbH



RENODRAIN NEPHROSTOMY SET



Set acc. Dr. Schüller

SKATER® EXCHANGE CATHETER

The set consist of:

- ✓ Pigtail-catheter made of polyurethane, distal part with hydrophilic coating

Characteristics of the catheter:

- ✓ exchange catheter for percutaneous drainage of abscesses, hematomas, ascites, seromas, nephrostomies, urinomas, bilomas, thoracentesises, cholecystotomies and other fluids
- ✓ soft polyurethane for kink-resistance and configuration memory
- ✓ hydrophilic coating for easy catheter insertion
- ✓ large oval drainage holes and large lumen are designed for maximal drainage capacity
- ✓ improved patient comfort due to minimal tissue trauma and reduced ulcerations

dimensions	
length	25,35 cm
diameter	6-14 CH
box	5 pcs

product code
08751106025
08751106035
08751107025
08751107035
08751108025
08751108035
08751110025
08751110035
08751112025
08751112035
08751114025
08751114035

CH cm

! "#\$ %& '()*+,-./:;<=>?@A#B CDEFGHIJKL MNOPQRSTUVWXYZ[\]^_`{|}~



SKATER® EXCHANGE CATHETER



Skater® exchange catheter

SKATER® ONE STEP PUNCTURE DRAINAGE SET

The set consist of

- ✓ trocar with Choice Lock™ snap lock
- ✓ Pigtail-catheter made of polyurethan

SKATER® ONE STEP PUNCTURE DRAINAGE SET, PIGTAIL CATHETER

Characteristics of the catheter:

- ✓ application during abscesses, ascites, thoracentesises, cholecastomies, hematomas, nephrostomies and other drainge procedures
- ✓ soft polyurethane for kink-resistance and configuration memory
- ✓ hydrophilic coating for easy catheter insertion
- ✓ large oval drainage holes and large lumen are designed for maximal drainage capacity
- ✓ improved patient comfort due to minimal tissue trauma and reduced ulcerations

dimensions catheter

length	20, 25 cm
diameter	6-14 CH
box	5 pcs

product code

08756006020
08756006025
08756007020
08756007025
08756008020
08756008025
08756010025
08756012025
08756014025

CH cm



SKATER® DRAINAGE SET



Skater® pigtail-catheter (1) with Choice Lock trocar (2+3)

URETERAL CATHETER STRAIGHT, CLOSED TIP

- ✓ radiopaque
- ✓ graduated in cm
- ✓ with stylet
- ✓ with injection adapter
- ✓ 2 lateral drainage holes

dimensions	
length	70 cm
diameter	3-7 CH
box	10 pcs

product code

08UK500370
08UK500470
08UK500570
08UK500670
08UK500770

CHcm

URETERAL CATHETER STRAIGHT, OPEN TIP

With same characteristics as the product above, additionally

- ✓ one side drainage hole
- ✓ adapted for 0,035" guidewire

product code

08UK600370
08UK600470
08UK600570
08UK600670
08UK600770

CHcm

URETERAL CATHETER TIEMANN, CLOSED TIP

With same characteristics as the product above, additionally

- ✓ 2 lateral drainage holes

product code

08UK510370
08UK510470
08UK510570
08UK510670
08UK510770

CHcm

URETERAL CATHETER TIEMANN, OPEN TIP

With same characteristics as the product above, additionally

- ✓ one side drainage hole
- ✓ adapted for 0,035" guidewire

product code

08UK540370
08UK540470
08UK540570
08UK540670
08UK540770

CHcm



URETERAL CATHETER



straight, closed
tip



straight, open
tip



Tiemann, closed
tip



Tiemann, open
tip

MIROMETER®

- ✓ closed urine measuring system for adults
- ✓ universal mount
- ✓ fixation clamp

dimensions	
tube length	120 cm
measuring chamber	500 ml
bag capacity	2000 ml
box	20 pcs.

product code
08M06050

MIROMETER® PÄD

- ✓ closed urine measuring system for children
- ✓ universal mount
- ✓ fixation clamp

dimensions	
tube length	120 cm
measuring chamber	200 ml
bag capacity	2000 ml
box	20 pcs.

product code
08M06030

MIROMETER® NEO 60

- ✓ closed urine measuring system for newborns
- ✓ with belt bracket

dimensions	
measuring chamber	60 ml
bag volume	150 ml
box	20 pcs.

product code
0520771



MIROMETER®



MiroMeter®



MiroMeter® Päd

MIROBAG® 2000
closed urine drainage system

- ✓ needle free urine sampling port
- ✓ universal mount
- ✓ fixation clamp and cord
- ✓ drain tube with stopcock. One hand operation for hygienic evacuation of the bag
- ✓ bacteria-proof ventilation filter at the catheter connector for fast urine flow into the bag reservoir
- ✓ drip chamber with integrated antireflux valve prevents the reflux of urine

dimensions	
tube length	120 cm
bag capacity	2000 ml
box	20 pcs.

product code
08M06032

MIROBAG® BIG 4000
closed urine drainage system for postoperative flushing

same properties as above and additional

- ✓ 4000 ml bag volume
- ✓ extra strong bed suspension

dimensions	
tube length	120 cm
bag capacity	4000 ml
box	20 pcs.

product code
08M06202

**MIROBAG®**

MiroBag®



MiroBag® Big

! #&! %&() *+,-./:;<=>?@# 71* #15 B0#%&0DD2 #2 %&8% EF! - "G 2 #2 %&8% H%#% N|>=JK@JLMLR4

Shaped Bandages for leg bag fixation.

- ✓ applicable with all leg bags up to 750 ml
- ✓ can be reused several times and can be washed (40° C)
- ✓ high wearing comfort
- ✓ no contact of the bag with skin
- ✓ prevents skin irritations
- ✓ reliable fixation

- ✓ individual adjustment in 4 sizes
- ✓ cushioning of fluid sounds

Dimensions	
length	S, M, L, XL
box	10 pcs

product code	
S	08745767001
M	08745767002
L	08745767003
XL	08745767004

COMFORTFIX LEG BAG PANT


- ✓ wide waist band
- ✓ skin-friendly and pollution-free organic tissue
- ✓ fits on both legs
- ✓ applicable with all common leg bag types
- ✓ easy to evacuate the bag
- ✓ can be washed up to 95°C
- ✓ available in 5 sizes

product code	
sizes	S, M, L, XL, XXL
box	1 pcs

product code	
S (60–67 cm)	087180
M (67–77 cm)	087180M
L (77–85 cm)	087181
XL (85–94 cm)	087182
XXL (95–105 cm)	087183



LEG BAG PANT ACCESSORIES



comfortfix leg bag pant

MiroValve® AUTOM. CATHETER VALVE FOR UROLOGY

- ✓ can't forget to open or close the valve
- ✓ automatically open once connected to the bag
- ✓ the valve closes automatically when bag is removed
- ✓ easy operation with one hand
- ✓ no dripping of urine
- ✓ applicable with all common catheters and bags
- ✓ permanent automatic indwelling by connection to the bag
- ✓ good for bladder training
- ✓ easy spray disinfection
- ✓ coloured operating elements
- ✓ valve system will be changed together with catheter
- ✓ can be used for biliary drainage too

All products and packaging are latex and pvc free.
A product video is provided here:
<http://www.youtube.com/watch?v=A88emF907pE>




MiroValve®

MiroValve® technical data

product code 089805

box 10 pcs



MiroValve® Automatic Catheter Valve

! # \$ % & ' () * + , - . / : ; < = > ? @ [\] ^ _ ` { | } ~ ¡ ¢ £ ¤ ¥ ¦ § ¨ © ª « ¬ ® ¯ ° ± ² ³ ´ µ ¶ · ¸ ¹ º » ¼ ½ ¾ ¿ À Á Â Ã Ä Å Æ Ç È É Ê Ë Ì Í Î Ï Ñ Ò Ó Ô Õ Ö × Ø Ù Ú Û Ü Ý Þ ß à á â ã

TOUHY-BORST ADAPTER

- ✓ with funnel connector

dimensions	
diameter	3-9 CH
box	10 pcs

product code	
	089935

TOUHY-BORST ADAPTER

for connecting an external drainage (Mono-J)

- ✓ female Luer Lock connector

dimensions	
diameter	3-9 CH
box	10 pcs

product code	
	089936

STEP CONE LUER-LOCK

- ✓ male¹
- ✓ female

ordering information	
box	10 pcs

product code	
male	089930
female	089931

UNIVERSAL CONNECTING TUBE

- ✓ silicone
- ✓ male Luer Lock connector

dimensions	
length	15, 30 cm
box	10 pcs

product code	
15 cm	089900
30 cm	089901

¹ is not illustrated



ACCESSORIES



Touhy-Borst with funnel



Touhy-Borst



step cone female



connecting tube

CATHETER EXCHANGE SET

For catheters with balloon capacity of 30 ml. The set is sterile and consists of:

- ✓ 1 kidney basin, pulp
- ✓ 1 plastic tray
- ✓ 5 gauze balls, plum-sized
- ✓ 1 gauze compress, 7,5 x 7,5 cm, eight-layer
- ✓ 1 drape with slit opening, 60 x 60 cm, two-layer
- ✓ 3 prefilled syringes with 10 ml aqua dest. glycerine solution (10 %)
- ✓ 1 empty syringe, 10 ml
- ✓ 1 fixing plaster Skinfix Maxi
- ✓ 1 pair of gloves, medium size, powder free
- ✓ 1 drape 60 x 60 cm, two-layer

ordering information

box	20 pcs
-----	--------

product code

08M5901

CATHETER STOPPER

ordering information

box	50 pcs
-----	--------

product code

089810

SKINFIX MAXI

- ✓ fixation plaster

ordering information

box	50 pcs
-----	--------

product code

089710



CATHETER EXCHANGE SET



Catheter exchange set



Catheter stopper